NADAL[®] RSV Test (test cassette)

REF 491008N-20

DE Gebrauchsanweisung	2	Symbols	11
EN Instructions for use	5	Our Teams	12

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Carl-Zeiss-Strasse 12
47445 Moers
Germany

Moers
Tel: +49 (2841) 99820-0
Fax: +49 (2841) 99820-1

Regensburg
Tel: +49 941 29010-0
Fax: +49 941 29010-50

www.nal-vonminden.com
info@nal-vonminden.com

Directors:
Sandra von Minden
Roland Meißner
Thomas Zander

Commercial reg. Kleve
HRB 5679
Steuer-Nr. 244/133/00130
UST-ID-Nr. DE 189 016 086

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Adenovirus

15. Referenzen

1. BARENFANGER et al., "Clinical and Financial Benefits of Rapid Detection of Respiratory Viruses: an Outcomes Study". Journal of Clinical Microbiology, August 2000, Vol 38 No 8, p. 2824-2828.

Rev. 0, 2021-10-29 OM

1. Intended Use

The NADAL® RSV Test is lateral flow chromatographic immunoassay for the qualitative detection of RSV antigens in human nasopharyngeal swab, wash or aspirate specimens. The test is intended for use as an aid in the diagnosis of RSV infections in patients showing typical symptoms of a respiratory tract infection. The test procedure is not automated and requires no special training or qualification. The NADAL® RSV Test is designed for professional use only.

2. Introduction and Clinical Significance

Although a wide variety of viral agents are capable of causing lower respiratory tract infections in children and adults, respiratory syncytial virus (RSV), influenza A and B viruses, parainfluenza viruses 1, 2, and 3, as well as adenovirus, are the most common. Of these, RSV and influenza A and B viruses are the most foremost causes of medically-attended acute respiratory illnesses. In addition to sharing a similar seasonal prevalence, it is important to remain aware that RSV as well as influenza A and B share overlapping clinical features and infection potential for certain high-risk patient groups (e.g. extremes of age, underlying cardiopulmonary disease or immunosuppression).

3. Test Principle

The NADAL® RSV Test enables the detection of RSV antigens through the visual interpretation of colour development on the internal test strip. Monoclonal anti-RSV antibodies are immobilised in the test line region (T) of the membrane. During the test, the specimen reacts with other anti-RSV antibodies which are conjugated to coloured particles and precoated onto the sample pad of the test cassette. The mixture then migrates along the membrane by capillary action and interacts with the reagents on the membrane. If there are a sufficient number of RSV antigens in the specimen, a red line will develop in the test line region (T) of the membrane. The presence of this line indicates a positive result, while its absence indicates a negative result. The formation of a green line in the control line region (C) serves as a procedural control, indicating that the proper volume of specimen has been added and membrane wicking has occurred.

4. Reagents and Materials Supplied

- 20 NADAL® RSV test cassettes, incl. disposable pipettes
 - 20 extraction tubes
- Additional material provided according to 93/42/EEC:
- 20 sterile swabs, CE0197

Jianguo Changfeng Medical Industry Co., Ltd
Touqiao Town, Guangling District, Yangzhou,
Jiangsu 225109 China (authorised EU
representative Llins Service & Consulting GmbH,
Obere Seegasse 34/2, 69124 Heidelberg,
Germany)

- 1 buffer (12 mL)
- 1 RSV positive control swab
- 1 package insert

5. Additional Materials Required

- Specimen collection containers (suitable for specimen material to be tested)

- Aspiration device with catheter and mucus trap (for aspirates)
- Timer
- Shaker or vortexer (if available)

6. Storage and Stability

Test kits should be stored at 2-30°C until the indicated expiry date. Test cassettes are stable until the expiry date printed on the foil pouches. Test cassettes must remain in the sealed foil pouches until use. Do not freeze test kits. Do not use tests beyond the expiry date indicated on the packaging. Care should be taken to protect test kit components from contamination. Do not use test kit components if there is evidence of microbial contamination or precipitation. Biological contamination of dispensing equipment, containers or reagents can lead to inaccurate results.

7. Warnings and Precautions

- For professional *in-vitro* diagnostic use only.
- Carefully read through the test procedure prior to testing.
- Do not use the test beyond the expiration date indicated on the packaging.
- Do not use test kit components if the primary packaging is damaged.
- Tests are for single use only.
- Do not add specimens to the reaction area (result area).
- In order to avoid contamination, do not touch the reaction area (result area).
- Avoid cross-contamination of specimens by using a new extraction tube for each specimen obtained.
- Do not substitute or mix components from different test kits.
- Do not use the buffer if it is discoloured or turbid. Discolouration or turbidity may be a sign of microbial contamination.
- Do not eat, drink or smoke in the area where specimens and test kits are handled.
- Wear protective clothing such as laboratory coats, disposable gloves and eye protection when specimens are being assayed.
- Handle all specimens as if they contain infectious agents. Observe established precautions for microbiological risks throughout all procedures and standard guidelines for the appropriate disposal of specimens.
- The test kit contains products of animal origin. Certified knowledge of the origin and/or sanitary state of the animals does not completely guarantee the absence of transmissible pathogenic agents. It is therefore recommended that these products be treated as potentially infectious and handled in accordance with usual safety precautions (e.g., do not ingest or inhale).
- Temperature can adversely affect test results.
- Used testing materials should be disposed of according to local regulations.

8. Specimen Collection and Preparation

Nasopharyngeal swab

- Bend the shaft of a sterile swab to follow the curve of the nasopharynx.

- Insert the swab into the nostril, pushing it gently into the posterior nasopharynx.
- Rotate the swab several times in order to obtain infected cells.
- For an optimal specimen, repeat the procedure in the other nostril.

Nasopharyngeal aspirate method (suction apparatus, sterile suction catheter):

- Add several drops of saline solution to each nostril.
- Insert a catheter through a nostril into the posterior nasopharynx.
- Apply gentle suction. Using a rotating motion, slowly withdraw the catheter.
- For an optimal specimen, repeat the procedure in the other nostril.


Note: Specimens should be tested as soon as possible after collection.

9. Test Procedure

Bring tests, specimens, buffer and/or controls to room temperature (15-30°C) prior to testing.


1a. Nasopharyngeal wash/aspirate

- I. Holding a pipette vertically, dispense 6 drops of the wash/aspirate specimen into the extraction tube.
- II. Holding the buffer bottle vertically, add 9 drops of buffer to the extraction tube.
- III. Mix the solution with a shaker, or vigorously by hand, for 1 minute.

1b Nasopharyngeal swab

- I. Holding the buffer bottle vertically, dispense 15 drops of buffer into the extraction tube
- II. Insert the swab with the collected nasopharyngeal specimen into the tube. Swirl the swab to mix and extract as much liquid as possible from it. Remove the swab from the tube.

2. Remove the test cassette from the foil pouch and use it as soon as possible. The best results will be obtained if the test is performed immediately after opening the foil pouch. Label the test cassette with the patient or control identification.

3. Place the test cassette on a clean and level surface.

4. Holding a pipette vertically, transfer exactly 4 drops of extracted specimen to the specimen well (S) of the test cassette.

Note: If the test does not run due to the specimen type, stir the added specimen in the specimen well (S)

with the pipette. If this does not work, add 1 or maximum 2 drops of buffer to the specimen well (S) until you observe a coloured liquid migrating along the membrane.

5. Start the timer.

6. Wait for the coloured line(s) to appear. Read the result after exactly 10 minutes. Do not interpret the result after more than 10 minutes.

10. Result Interpretation

Positive:

A green line develops in the control line region (C) and another red line develops in the test line region (T).

Note: The colour intensity in the test line region (T) may vary depending on the analyte concentration present in the specimen. Any shade of colour in the test line region (T) should be considered positive. Note that this is a qualitative test only and it cannot determine the analyte concentration in the specimen.

Negative:

A green line develops in the control line region (C). No red line develops in the test line region (T).

Invalid:

The green control line (C) fails to develop. Results from any test which has not produced a control line at the specified reading time must be discarded.

Please review the procedure and repeat the test with a new test cassette. If the problem persists, discontinue using the test kit immediately and contact your distributor.

Insufficient specimen volume, incorrect operating procedure or expired tests are the most likely reasons for control line failure.

11. Quality Control

The internal procedural control is included in the test cassette:

A green line appearing in the control line region (C) is considered an internal positive procedural control. It confirms sufficient specimen volume, correct procedural technique and adequate membrane wicking.

Good laboratory practice (GLP) recommends the use of control materials to ensure proper test kit performance.

An RSV positive control swab is provided with each test kit to monitor test performance in the detection of RSV antigens. To perform control testing, follow the same procedure as with a nasopharyngeal swab specimen.

12. Limitations

- The NADAL® RSV Test is for professional *in-vitro* diagnostic use only.
- The NADAL® RSV Test only indicates the presence of RSV antigens in the specimen and should be used for the

qualitative detection of RSV antigens in nasopharyngeal swabs, washes and aspirates only. Neither the quantitative value nor the rate of increase in RSV antigen concentration can be determined with this test.

- If the test result is negative but clinical symptoms persist, additional testing using other clinical methods is recommended. A negative result does not at any time preclude the possibility of RSV infection.
- The test provides a presumptive diagnosis of RSV infection. As with all diagnostic tests, all results should be interpreted in conjunction with other clinical information and laboratory findings available to the physician.

13. Expected Values

RSV is generally considered the most frequent cause of pneumonia, bronchiolitis and tracheobronchitis in infants and young children. It is known to be the etiologic cause in 14-27% of cases of pneumonia in the elderly during the winter season.

14. Performance Characteristics

Sensitivity and specificity

Different virus extract dilutions were tested directly in buffer or added to negative nasopharyngeal specimens.

The results yield a sensitivity of 95% and specificity of >99% in comparison with another commercially available rapid test Alere BinaxNOW® RSV.

Cross-reactivity

An evaluation was performed to determine the cross-reactivity of the NADAL® RSV Test. There was no cross-reactivity found with common respiratory pathogens, other organisms or substances occasionally present in nasopharyngeal specimens:

Influenza type A virus

Influenza type B virus



Adenovirus

15. References

1. BARENFANGER et al., "Clinical and Financial Benefits of Rapid Detection of Respiratory Viruses: an Outcomes Study". Journal of Clinical Microbiology, August 2000, Vol 38 No 8, p. 2824-2828.

Rev. 0, 2021-10-29 OM

Symbol	Deutsch	English	Français	Español	Italiano	Polski
	CE Konformitätszeichen	CE marking of conformity	Conforme aux normes européennes	Conformidad europea	Conformità europea	Znak zgodności CE
	Gebrauchsanweisung beachten	Consult instructions for use	Consulter la notice d'utilisation	Consúltense las instrucciones de uso	Consultare le istruzioni per l'uso	Przestrzegać instrukcji obsługi
	<i>in-vitro</i> -Diagnostika	<i>in-vitro</i> diagnostic medical device	Dispositif médical de diagnostic <i>in-vitro</i>	Producto sanitario para diagnóstico <i>in-vitro</i>	Dispositivo medico-diagnostico <i>in-vitro</i>	Tylko do diagnostyki <i>in-vitro</i>
	Temperaturbegrenzung	Temperature limitation	Limites de température	Limite de temperatura	Limiti di temperatura	Temperatura przechowywania
	Chargenbezeichnung	Batch code	Numéro de lot	Código de lote	Codice lotto	Numer serii
	Nicht zur Wiederverwendung	Do not reuse	Ne pas réutiliser	No reutilizar	Non riutilizzare	Tylko do jednorazowego użytku
	Verwendbar bis	Use by	Utiliser jusqu'au	Fecha de caducidad	Utilizzare entro	Data ważności
	Bestellnummer	Catalogue Number	Référence du catalogue	Número de catálogo	Riferimento di Catalogo	Numer katalogowy
	Hersteller	Manufacturer	Fabricant	Fabricante	Fabbricante	Producent
	Ausreichend für <n> Ansätze	Sufficient for <n> tests	Suffisant pour "n" tests	Suficiente para <n> utilizations	Sufficiente per "n" saggi	Wystarczający na <n> Powtórzeń

Symbol	Português	Český	Suomi	Svenskt	Nederlands	Dansk	Norsk
	Conformidade com as normas europeias	CE certifikát	CE-merkitty	CE-märkning	CE-markering	CE-mærkning	CE standardisert
	Consultar as instruções de utilização	Viz návod k použití	Katso käyttöohjetta	Läs bruksanvisningen	Raadpleeg de gebruiksaanwijzing	Se brugsanvisningen	Les bruksanvisning nøye
	Dispositivo médico para diagnóstico <i>in-vitro</i>	Diagnostický zdravotnický prostředek <i>in-vitro</i>	<i>in-vitro</i> -diagnostikkaan tarkoitettu lääkinnällinen laite	Medicinteknisk produkt avsedd för <i>in-vitro</i> -diagnostik	Medisch hulpmiddel voor <i>in-vitro</i> diagnostiek	Medicinsk udstyr til <i>in-vitro</i> -diagnostik	<i>in-vitro</i> diagnostic medisinsk enhet
	Limites de temperatura	Teplotní omezení	Lämpötilarajat	Temperaturbegränsning	Temperatuurlimiet	Temperaturbegrænsning	Temperaturbegrensning
	Código do lote	Kód šarže	Eräkoodi	Satsnummer	Code van de partij	Batchkode	Merking
	Não reutilizar	Pro jednorázové použití	Kertakäyttöinen	Får inte återanvändas	Niet opnieuw gebruiken	Må ikke genbruges	Må ikke brukes om igjen
	Prazo de validade	Spotřebuje do	Käytettävä viimeistään	Används före	Houdbaar tot	Udløbsdato	Tidtaking
	Número de catálogo	Katalogov číslo	Luettelonumero	Listnummer	Catalogus nummer	Bestil lingsnummer	Katalog nummer
	Fabricante	Výrobce	Valmistaja	Tilverkare	Fabrikant	Fabrikant	Produsent
	Suficiente para <n> test	Dostačuje pro <n> testů	Lukumäärä <n> test	Räcker till <n> test	Voldoende voor <n> test	Tilstrækkeligt til <n> test	Tilstrækkelig for <n> tester

Our Teams

Germany:

Regensburg

Tel: +49 941 290 10-0
Fax: +49 941 290 10-50

Moers

Tel: +49 2841 99820-0
Fax: +49 2841 99820-1

Austria:

Tel: +49 941 290 10-29
Free Tel: 0800 291 565
Fax: +49 290 10-50
Free Fax: 0800 298 197

UK & Ireland:

Tel: +49 941 290 10-18
Free Tel –UK: 0808 234 1237
Free Tel – IRE: 1800 555 080
Fax: +49 290 10-50

France:

France Tel: 0800 915 240
France Fax: 0800 909 493

Switzerland

Swiss Tel: 0800 564 720
Swiss Fax: 0800 837 476

Belgium

Belgium Tel: 0800 718 82
Belgium Fax: 0800 747 07

Luxembourg

Lux. Tel: 800 211 16
Lux. Fax: 800 261 79

Spain:

Tel: +49 941 290 10-759
Free Tel: 900 938 315
Fax: +49 941 290 10-50
Free Fax: 900 984 992

Italy:

Tel: +49 941 290 10-34
Fax: +49 941 290 10-50

Poland:

Tel: +49 941 290 10-44
Free Tel: 00 800 491 15 95
Fax: +49 941 290 10-50
Free Fax: 00 800 491 15 94

Portugal:

Tel: +49 941 290 10-735
Tel. Verde: 800 849 230
Fax: +49 941 290 10-50
Fax Verde: 800 849 229

Netherlands:

Tel: +31 30 75 600
Free Tel: 0800 0222 890
Fax: +31 70 30 30 775
Free Fax: 0800 024 9519

Nordic countries (Finland, Norway, Sweden, Denmark):

Tel: +31 703075 607
Free Tel: +45 80 88 87 53
Tax: +31 703030 775

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Carl-Zeiss-Strasse 12 • 47445 Moers • Germany

www.nal-vonminden.com • info@nal-vonminden.com

Tel: +49 2841 99820-0 • Fax: +49 2841 99820-1

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nal von minden GmbH